


FLOCARE INFINITY PUMP SHORT MANUAL

Preparing and setting up the pump



Open the pump door by pressing gently up on the lower wall (A) and rotating the door upwards (B) at the same time.

 Position the looped section of the silicone tubing around the rotor. Stretch slightly to seat the cassette into the pump.

Close the pump door

 you should hear a
 small click when it is
 closed properly. The
 pump is now ready to
 be switched on.

Operating instructions

Clearing pump memory and priming the feeding set

ON/OFF

When switching the pump on or off, press this key for 2 seconds.

INFO

CLR

This will clear the previous volume.
The rate and dose will remain unchanged.

FILL SET

Pressing this key for 2 seconds will automatically fill the feeding set. Connect the feeding set to the feeding tube and choose appropriate feeding program 1 or 2.

- 1 Continuous feeding
- ml/h
- + 400ml/h
- DOSE=VOL
- CLR
- START/STOP

- Pre-set dose
- ml/h
- + 400ml/h
- DOSE=VOL





START/STOP

PROBLEM SOLVING INSTRUCTIONS

If your pump alarms:

STEP 1 Look at the pump and note the error message displayed on the screen.

STEP 2 Push the **START/STOP** button to stop the alarm noise.

STEP 3 Refer to the chart below for steps to resolve the problems.



Error messages and instructions

Please refer to the detailed **Instructions for Use** booklet that comes with your pump for more information.

Problem	Instructions
NO SET	Press ON/OFF to switch off the pump. Insert the Flocare® Infinity™ set into the pump following instructions for use. Press ON/OFF to restart the pump.
PUSH START	The pump has been untouched in hold mode for over 3 minutes. Extend the hold mode by a further 3 minutes by pressing START/STOP . Check the pump programming is correct the pump and press START/STOP to start the pump.
END OF DOSE	The programmed dose/volume has been administered. Switch off the pump by pressing ON/OFF.
PROG	This could be a result of holding your finger down on the CLR button for too long, or no flow rate has been installed. Make sure the correct flow rate is programmed.
DOOR	The door is not correctly closed. Make sure the door is properly closed prior to starting a feeding program.
AIR	Press START/STOP . Check if the set is inserted correctly and check the air sensor. Remove the air in the set by pressing the FILL SET key once. Once all visible air has been removed from the feeding set, press FILL SET to pause the pump movement. Recommence feeding by pressing START/STOP .
BATT	Connect the charging cable to the pump, and plug the adaptor in to an electrical outlet. Ensure that the adaptor (the part that connects to the electrical outlet) is clicked in to place, sitting flat and flush - a green light should display when properly connected. A full charge takes approximately 6 hours. The pump can be used whilst on charge.
BATT E + F	Battery failure. Please contact the pumps team to organise a pump service (contact information below).
OEE IN	Press START/STOP to silence the alarm. Check the feeding set for any kinks in the tubing or visible blockages. If required, remove the set from the pump and clean the pump sensors before reinserting the cartridge back in to place. Press START/STOP to recommence feeding.
LOEK	The programming feature is locked in the set-up mode of the pump. Refer to instruction booklet to unlock.
ERO) - ER99	Electronic error. Switch the pump off, make sure the pump door is closed and turn the pump back on. If the error persists, please refer the pump for service by contacting Nutricia.
NO PLUG SYMBOL	First check the adaptor is firmly pressed down on the area housing the two pins, you should hear a click when this is in place. Reconnect to the wall. If there is still no plug symbol on the pump and no green light on the adaptor, try using another wall outlet. If the pump is still not able to charge, request a replacement of the adaptor.

Nutricia Pumps Team - 1800 060 051 (Select option 4 for the pumps team) Mon-Fri: 9am-5pm AEST[^]

