

FLOCARE® ENTERAL FEEDING PUMPS

Flocare® Infinity™ II and Flocare® Infinity™+ (Plus) Pumps are small, portable enteral feeding pumps offering maximum flexibility and mobility. Flocare® Enteral Feeding Pumps are designed for controlling the delivery rate of enteral feeds into the gastrointestinal tract.

FEATURES

Flocare® Infinity™ II Pump

- One-way cassette that uses rotary peristaltic technology to prevent free-flow of feed.
- Closed pack system.
- Automatic priming.
- Extensive alarm options (upstream and downstream occlusion detection, air bubble detection, battery capacity, end of dose).
- Lightweight (392g).
- Size: 144 x 105 x 50mm.
- Built in continuous dual memory of feeding program and volume delivered for easy programming.
- Clear LED display.
- Flow rate 1-400ml/hr (in 1 ml increments).
- Dose feature 1-4000ml (in 1 ml increments).
- Lithium ion battery 3.6V DC 2000 mAh, 6-hour charging time.
- 24-hour battery life on full charge, at 125ml/hr.
- Compatible with Flocare® Infinity Feeding Sets.
- Easy to clean, water resistant.
- Access to hidden menu with: ability to lock feeding regiment to prevent feeding rate or dose being altered; alarm on/off options; alarm volume options and LCD light on/off option.

Flocare® Infinity™+ (Plus) Pump

- All Flocare® Infinity™ II Pump features included.
- Additional function to enable bolus feeding at timed intervals.
 - Volume: 1-500ml (in 1 ml increments).
 - Time intervals: 0.5-12 hours (in 30 min intervals).
 - Flow rate: 400ml/hr or 600ml/hr.

Flocare® Infinity™ II Pump



Product Code: 40407

Units Per Carton: 1

Flocare® Infinity™+ (Plus) Pump



Product Code: 40408

Units Per Carton: 1

Directions for Use

- Service check is required at least every 2 years.
- See Nutricia's 'A step-by-step guide to setting up Nutricia's ready to hang system' for detailed directions for use.

Storage

- Store in a clean, dry environment.
- Protect from dust, rain and direct sunlight.

Ordering Information

To order contact Nutricia Customer Care **1800 889 480**.