




FLOCARE® ACCESSORIES

The Flocare® accessories range includes feeding connectors, adaptors, containers, and pump accessories to aid enteral feeding. All Flocare® Feeding Tubes, Flocare® Feeding Sets and ENFit™ enteral syringes are ENFit™ compatible.

Adaptors	Product description	Code	Units per carton
<p>ENFit™ Transition Adaptor</p> 	<p>Designed to be inserted into a non-ENFit™ feeding tube to enable connection to an ENFit™ feeding set.</p> <p>Recommended for single patient use for the lifespan of the feeding tube.</p>	122073	50
<p>Flocare® Transition Step Connector</p> 	<p>Designed to be attached to an ENFit™ feeding set to enable connection to a non-ENFit™ feeding tube.</p> <p>Recommended for 24-hour, single use, single patient only.</p>	89828	30

Containers	Product description	Code	Units per carton
<p>Flocare® Container 0.5L</p> 	<p>Designed to be used to deliver decanted product, reconstituted powder based formula and/or water via an enteral feeding pump or gravity feeding set.</p> <p>Recommended for 24 hour, single use only.</p>	40441	10
<p>Flocare® Container 1L</p> 	<p>Designed to be used to deliver decanted product, reconstituted powder based formula and/or water via an enteral feeding pump or gravity feeding set.</p> <p>Recommended for 24 hour, single use only.</p>	40455	10
<p>Flocare® Top-fill Container 1.3L</p> 	<p>Designed to be used to deliver decanted product, reconstituted powder based formula and/or water via an enteral feeding pump or gravity feeding set.</p> <p>Top-fill design for ease of use.</p> <p>Recommended for 24 hour, single use only.</p>	70139	10

To order, please call:
Nutricia Customer Care 1800 889 480

Mobile Accesories

Flocare® Universal Frame



Compatible with:

- Flocare® Infinity II and Flocare® Infinity III Pumps
- Full range of Nutricia OpTri Bottles and Flocare® Containers
- Giving Sets – with and without drip chamber
- Existing Flocare® backpacks
 - Suitable for stationary and mobile usage
 - Adjustable shelf to accommodate different sized bottles and containers.

Code

188229

Units per carton

1

Flocare® Infinity GO Bag



Designed to optimise patient mobility and product protection when tube feeding.

Capacity to carry both the Flocare® Infinity Pump, enteral feed and Universal Frame at the same time.

Code

44497

Units per carton

1

For more information, please call:
Nutricia Clinical Care Line 1800 060 051

NUTRICIA
LIFE-TRANSFORMING NUTRITION

Nutricia Australia Pty Ltd
The Zenith, Tower B, Level 12,
821 Pacific Hwy, Chatswood NSW 2067
www.nutricia.com.au/adult

NC4656-Nov23 / NUT1517