




# FLOCARE® ACCESSORIES




The Flocare® accessories range includes feeding connectors, adaptors, containers, and pump accessories to aid enteral feeding.

All Flocare® Feeding Tubes, Flocare® Feeding Sets and ENFit™ enteral syringes are ENFit™ compatible.

ADAPTORS	PRODUCT DESCRIPTION	CODE	UNITS PER CARTON
 <p data-bbox="225 1055 536 1088"><b>ENFit™ Transition Adaptor</b></p>	<p data-bbox="644 831 1023 1025">Designed to be inserted into a non-ENFit™ feeding tube to enable connection to an ENFit™ feeding set. Recommended for single patient use for the lifespan of the feeding tube.</p>	<p data-bbox="1091 913 1182 943">122073</p>	<p data-bbox="1310 913 1347 943">50</p>
 <p data-bbox="181 1435 576 1469"><b>Flocare® Transition Step Connector</b></p>	<p data-bbox="671 1211 995 1406">Designed to be attached to an ENFit™ feeding set to enable connection to a non-ENFit™ feeding tube. Recommended for 24-hour, single use, single patient only.</p>	<p data-bbox="1094 1294 1182 1323">89828</p>	<p data-bbox="1310 1294 1347 1323">30</p>
 <p data-bbox="197 1816 563 1850"><b>Flocare® ENFit™ Bolus Adaptor</b></p>	<p data-bbox="644 1592 1023 1787">Designed to be inserted into ENPlus enteral feeding packs to enable connection to an ENFit™ syringe to allow for bolus feeding. Recommended for 24-hour, single tube feed pack use.</p>	<p data-bbox="1094 1675 1182 1704">89740</p>	<p data-bbox="1310 1675 1347 1704">30</p>

To order, please call:

**Nutricia Customer Care 1800 889 480**

CONTAINERS	PRODUCT DESCRIPTION	CODE	UNITS PER CARTON
 <p>Flocare® Container 0.5L</p>	Designed to be used to deliver decanted product, reconstituted powder based formula and /or water via an enteral feeding pump or gravity feeding set. Recommended for 24 hour, single use only.	210440	10
 <p>Flocare® Container 1L</p>	Designed to be used to deliver decanted product, reconstituted powder based formula and /or water via an enteral feeding pump or gravity feeding set. Recommended for 24 hour, single use only.	210577	10
 <p>Flocare® Top-fill Container 1.3L</p>	Designed to be used to deliver decanted product, reconstituted powder based formula and /or water via an enteral feeding pump or gravity feeding set. Top-fill design for ease of use. Recommended for 24 hour, single use only.	70139	10

MOBILE ACCESSORIES	PRODUCT DESCRIPTION	CODE	UNITS PER CARTON
 <p>Flocare® Infinity GO Frame</p>	Table-top, lightweight, plastic frame that easily fixes to the Flocare® Infinity pump and will allow enteral pack / bottle, Op Tri bottles or Flocare® containers to be attached, for ease of delivery.	43247	1
 <p>Flocare® Infinity GO Bag</p>	Designed to optimise patient mobility and product protection when tube feeding. Capacity to carry both the Flocare® Infinity pump, enteral feed and GO frame at the same time.	44497	1

For more information, please call:  
**Nutricia Clinical Care Line 1800 060 051**

**NUTRICIA**  
LIFE-TRANSFORMING NUTRITION

Nutricia Australia Pty Ltd.  
Lvl 4, Building D, 12-24 Talavera Road,  
Macquarie Park, NSW, 2113.  
[www.nutriciamedical.com.au](http://www.nutriciamedical.com.au)

NC4656-Apr20 / NUT1141